



# Gorski Hip Cup Extraction Hook

Designed by Jerrold Gorski, MD

## Helps in the removal of a hip cup

Designed to quickly fit into a screw hole of a hip cup after the screws have been removed and the cup loosened. The slap hammer helps to remove the cup in the angle it was inserted.

| PRODUCT NO'S:                     |   |
|-----------------------------------|---|
| <b>Hook for 6.5mm Screw Holes</b> |   |
| 3660                              | [Hook with Standard (#3925) Slap Hammer]        |
| 3660-01                           | [Hook w/o Slap Hammer]                          |
| <b>Hook for 5.0mm Screw Holes</b> |   |
| 3665                              | [Hook with Standard (#3925) Slap Hammer]        |
| 3665-01                           | [Hook w/o Slap Hammer]                          |
| <b>Optional Items:</b>            |   |
| 3935                              | [Extra Large Slap Hammer] Thread Gauge: 3/8"-16 |



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

1.800.548.2362