

Gupta Disc Space Spreaders with Easy Release Locking Mechanism

Designed by Munish C. Gupta, MD


Designed to distract open collapsed disc spaces

Locking ratchet mechanism helps prevent accidental release, and provides for controlled adjustment and easy release.



Standard Straight

Large Straight



Standard Curved

Large Curved

PRODUCT NO'S:

4290 [Standard Straight]
Overall Length: 8.5" (21,6 cm)
Blade Width: 10 mm
Blade Thickness - Closed: 3 mm
Opens to: 22 mm

4292 [Standard Curved]
Overall Length: 9.25" (23,5 cm)
Blade Width: 10 mm
Blade Thickness - Closed: 3 mm
Opens to: 22 mm



4291 [Large Straight]
Overall Length: 11.5" (29,2 cm)
Blade Width: 13 mm
Blade Thickness - Closed: 4 mm
Opens to: 25 mm

4293 [Large Curved]
Overall Length: 12.5" (31,8 cm)
Blade Width: 13 mm
Blade Thickness - Closed: 4 mm
Opens to: 25 mm

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362