

Gupta Extended Osteotome

Designed by Munish C. Gupta, MD

Designed to help cut bone and cartilage in procedures such as facetectomies and vertebrectomies



PRODUCT NO:

5233

Overall Length: 11" (27,9 cm)
Blade Width: .5" (13 mm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362