


Gupta Probe Set

Designed by Munish C. Gupta, MD

A universal spine probe set designed to help cannulate pedicles, with various sized and shaped tips for use in cervical, thoracic, lumbar, and sacral pedicles

PRODUCT NO'S:	
5005-00	[Gupta Probe Set]
Set Includes / Available Individually:	
5005-01	[Spine Probe, Straight] Overall Length: 9.43" (24 cm) Handle Length: 5" (12,7 cm) Probe Shaft Diameter Tapers Down To: 1 mm Probe Ball End Diameter: 1.5 mm Depth Markings at: 2, 3, 4, and 5 cm
5005-02	[Spine Probe, Bent] Overall Length: 9.1" (23,1 cm) Handle Length: 5" (12,7 cm) Probe Shaft Diameter Tapers Down To: 1.5 mm Probe Ball End Diameter: 2.8 mm Depth Markings at: 2, 3, 4, and 5 cm
5005-03	[T-Handle Spine Probe, Flat Tip] Overall Length: 10.25" (26 cm) T-Handle Width: 2.75" (7 cm) Handle Length: 5" (12,7 cm) Probe Shaft Diameter: 3.6 mm Flat Tip Length: 20 mm Flat Tip Width: 1.2 mm Depth Markings at: 2, 3, 4, 5, and 6 cm
5005-04	[Pointed Thoracic Pedicle Probe] Overall Length: 9.43" (24 cm) Ball Handle End Diameter: 1.5" (3,8 cm) Handle Length (including ball end): 6.4" (16,3 cm) Before Tip, Probe Shaft Diameter Tapers Down To: 3 mm Depth Markings at: 1, 2, 3, 4, and 5 cm
5005-05	[Rounded Lumbar Pedicle Probe] Overall Length: 9.09" (23 cm) Ball Handle End Diameter: 1.5" (3,8 cm) Handle Length (including ball end): 6" (15,3 cm) Probe Shaft Diameter Tapers Down To: 4 mm Depth Markings at: 3, 4, 5, and 6 cm
5005-06	[Pediatric Pedicle Probe] Overall Length: 8.45" (21,5 cm) Ball Handle End Diameter: 1.5" (3,8 cm) Handle Length (including ball end): 5" (12,7 cm) Before Tip, Probe Shaft Diameter Step Tapers From: 5 to 2 mm Depth Markings at: 1, 2, 3, 4, 5, and 6 cm



Spine Probe, Straight
#5005-01

Spine Probe, Bent
#5005-02

Handle Spine Probe, Flat Tip
#5005-03

Pointed Thoracic Pedicle Probe
#5005-04

Rounded Lumbar Pedicle Probe
#5005-05

Pediatric Pedicle Probe
#5005-06

www.innomed.net



INNOMED

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.