


# Harvey Lumbar Bone Graft Sled Assembly

Designed by Charles Harvey, DO

*Designed to help deliver and tamp morselized bone graft to transverse processes during lumbar spinal fusion*



**PRODUCT NO'S:**

5083-00 [Harvey Lumbar Bone Graft Assembly]

**Set Includes / Available Individually:**

5083-01 [Harvey Lumbar Bone Graft Sled]

Overall Length: 10.625" (27 cm)

Sled Width: 1.5" (3,8 cm)

5083-02 [Harvey Lumbar Bone Graft Pusher]

Overall Length: 11.5" (29,2 cm)

Pusher End Width: 1.3" (3,3 cm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe**

Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**

**INNOMED** <sup>®</sup>