

Hibbs Retractors

Designed for soft tissue retraction by either the toothed end or curved handle end



The standard and medium sizes are useful in anterior approach THA to retract the Tensor Fascia Lata and Rectus Femoris during the initial approach to the anterior hip capsule. Helps expose the Lateral Circumflex vessels.

The extra large size can be used in large patients when more leverage and depth is needed.

PRODUCT NO'S:
6230 [Extra Large] Overall Length: 14.25" (36,2 cm) Handle Length: 13" (33 cm) Blade Depth: 6.5" (16,5 cm) Blade Width: 38 mm
6235 [Medium] Overall Length: 10.75" (27,3 cm) Handle Length: 9.75" (24,8 cm) Blade Depth: 4.5" (11,4 cm) Blade Width: 25 mm
6240 [Standard] Overall Length: 8.75" (22,5 cm) Handle Length: 8" (20,3 cm) Blade Depth: 3" (7,6 cm) Blade Width: 25 mm



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362