

Hibbs Retractors

Designed for soft tissue retraction by either the toothed end or curved handle end

Extra large used in large patients when more leverage and depth is needed.

Product No's:

6230 [Extra Large]
Overall Length: 14.25"
Handle Length: 13"
Blade Depth: 6.5"
Blade Width: 38mm

6235 [Medium]
Overall Length: 10.75"
Handle Length: 9.75"
Blade Depth: 4.5"
Blade Width: 25mm

6240 [Standard]
Overall Length: 8.75"
Handle Length: 8"
Blade Depth: 3"
Blade Width: 25mm



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

