


# Humeral Protection Plates

Designed by Ronald E. Delanois, MD

*Helps protect the proximal humerus from fracture after humeral head osteotomy*

Plate is placed on the proximal humerus after the initial osteotomy of the humeral head for total shoulder replacement. Helps protect the proximal humerus from fracture as the humerus is retracted to gain visualization of the glenoid to prepare it for a glenoid implant.



**PRODUCT NO'S:**

5259-01 [46mm]

5259-02 [50mm]

**PROUDLY  
MADE  
IN THE  
USA**

ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**INNOMED**

**1.800.548.2362**