


Iliac Screw Retractor Blade

Designed by Munish C. Gupta, MD

Designed to help assemble a long spinal construct to the pelvis, and can also help anyone performing a lumbosacral construct to the pelvis



PRODUCT NO:
7450-08C
Overall Length: 7.625" (19,4 cm)
Handle Length: 6.5" (16,5 cm)
Blade Depth: 4.35" (11 cm)
Blade Width: 1.75" 45 mm
Prong Width at Tip: 4 mm



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2025 Innomed, Inc.