






Incavo Wire Passer

Designed by Stephen J. Incavo, MD

Designed to pass multiple cerclage wires around a bone during a multiple wire wrap procedure



PRODUCT NO'S:
8610-01 [Small] Overall Length: 7.5" (19,1 cm) Accepts Wire Up To: 4 mm (5/32") For Bone Diameter Up To: 1.2" (3 cm)
8610-02 [Large] Overall Length: 8.675" (21,9 cm) Accepts Wire Up To: 4 mm (5/32") For Bone Diameter Up To: 2.4" (6 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362