

Inferior Acetabular Retractors

Help provide better access to the intramedullary canal



Standard



Deep Standard



Narrow

PRODUCT NO'S:
6250 [Standard] Overall Length: 12" (30,5 cm) Handle Length: 8" (20,3 cm) Blade Height Above Prongs: 3" (7,6 cm) Blade Width: 51 mm Prong Width: 5,1 mm 9,7 mm Gap 5,1 mm
6250-02 [Deep Standard] Overall Length: 12,5" (31,8 cm) Handle Length: 8" (20,3 cm) Blade Height Above Prongs: 4" (10,2 cm) Blade Width: 51 mm Prong Width: 5,1 mm 9,7 mm Gap 5,1 mm
6255 [Narrow] Overall Length: 12" (30,5 cm) Handle Length: 8" (20,3 cm) Blade Height Above Prongs: 3,25" (8,3 cm) Blade Width: 32 mm Prong Width: 5,1 mm 9,7 mm Gap 5,1 mm



Scan to
Launch Our
Website



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71



1.800.548.2362