

# Infero-Posterior Acetabular Retractor with Modular Handle – Left and Right

*Designed to be placed with the point at 6 o' clock and the retractor's axilla resting on the ischium, while the wing blade is used to retract the remaining capsule from the posterior lip of the acetabulum, and the optional screw-in modular handle can be used for additional leverage and maneuverability*



<b>PRODUCT NO'S:</b>
C1007-H-00 [Infero-Posterior Acetabular Retractor Set – Left]
<b>Set Includes/ Available Separately:</b>
C1007 [Infero-Posterior Acetabular Retractor Only – Left] Overall Length: 14" (35,6 cm) Depth from Bend: 4.5" (11,4 cm) Fixed Handle Width: 5.5" (14 cm)
C1006 [Modular Handle] Overall Length: 4.875" (12,4 cm) Handle Length: 4.5" (11,4 cm)



<b>PRODUCT NO'S:</b>
C1008-H-00 [Infero-Posterior Acetabular Retractor Set – Right]
<b>Set Includes/ Available Separately:</b>
C1008 [Infero-Posterior Acetabular Retractor Only – Right] Overall Length: 14" (35,6 cm) Depth from Bend: 4.5" (11,4 cm) Fixed Handle Width: 5.5" (14 cm)
C1006 [Modular Handle] Overall Length: 4.875" (12,4 cm) Handle Length: 4.5" (11,4 cm)



ISO 13485:2016  
**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**