


Johnson Low Profile Foot & Ankle Retractors

Designed by Michael Johnson, MD

Designed for soft tissue retraction in the foot and ankle



PRODUCT NO'S:	
1636-01 [Straight Handle]	Overall Length: 7.5" (19 cm) Blade Width: .41" (1 cm)
1636-02 [Double Bent Handle]	Overall Length: 8.875" (22,5 cm) Blade Width: .41" (1 cm)



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362