

# Kawell Short Army Navy Retractor

Designed by Ron Kane, DPM

*A short (4.75") handled Army Navy retractor,  
especially useful with a gastrocnemius recession*



**PRODUCT NO:**

**1148**

Overall Length: 4.75" (12,1 cm)  
Large End Blade Length: 1.75" (4,4 cm)  
Large End Blade Width: .625" (1,6 cm)  
Small End Blade Length: 1" (2,5 cm)  
Small End Blade Width: .3125" (0,8 cm)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

[www.innomed.net](http://www.innomed.net)



**INNOMED**



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**