


# Kenerly Double Parallel Femoral Neck Disc Grasper

Design modification by J. Lex Kenerly, III, MD of an original design by Alexander D. Rosenstein, MD

*Designed to remove the central disc of a double, parallel cut femoral neck osteotomy when performing THA*



**PRODUCT NO:**



**1720-01**

Overall Length: 10" (25,4 cm)

Jaw Dimensions: 1,44" x .72" (36,6 x 18,3 mm)

Lower Jaw Thickness: 1 mm

MADE EXCLUSIVELY  
FOR INNOMED IN  
GERMANY



[www.innomed.net](http://www.innomed.net)



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**