


Kirschenbaum Foot Positioners

Designed by Ira Kirschenbaum, MD

Helps eliminate the use of sand bags under the drape during total knee surgery

The foot rest is dome shaped for optimal foot contact and positioning the leg in flexion, and can be rotated. The unit can be used under the drape by attaching it to a standard table attachment or it can be sterilized for use on top of the drape. It can be attached to the table with the optional sterilizable table clamp. The foot rest is supplied with a removable, sterilizable silicone foot pad.



PRODUCT NO'S:
2590 [Long] Dimensions: 15.5" x 6" (39,4 cm x 15,2 cm)
2591 [Short] Dimensions: 9.5" x 6" (24,1 cm x 15,2 cm)
Optional & Replacement Parts:
2590-P [Long Replacement Pad] Dimensions: 16" x 9" (40,7 cm x 22,9 cm)
2591-P [Short Replacement Pad] Dimensions: 9.5" x 9.25" (24,1 cm x 23,5 cm)
2595 [Optional Table Clamp]



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362