


# Knee Retractors with Easy Grip Handles

*Helps provide excellent visibility and  
ligament protection during total and  
unicondylar knee replacement surgery*

Silicone handles help reduce holding fatigue.



Small Hohmann Retractor



Condylar Retractor



Soft Tissue Retractor



Superior Retractor

| PRODUCT NO'S:  |
|--|
| SS3035 [Small Hohmann Retractor]<br>Overall Length: 7" (17,8cm)<br>Blade Width: 25mm |
| SS3037 [Condylar Retractor]<br>Overall Length: 7" (17,8cm)<br>Blade Width: 12mm      |
| SS3038 [Superior Retractor]<br>Overall Length: 8.25" (21cm)<br>Blade Width: 31mm     |
| SS3042 [Soft Tissue Retractor]<br>Overall Length: 8.25" (21cm)<br>Blade Width: 36mm  |



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**