

# Kodkani Tissue Elevator Suture/Graft Passer

Designed by Pranjai Kodkani, MD

*Designed for MPFL reconstruction basket weave technique, and helpful for mini-open ligament reconstruction surgeries for graft passage*

Can also be used for:

- ▶ Periosteum/soft tissue elevator or freer
- ▶ Percutaneous passage of tendon/ligament graft/suture
- ▶ Stripping tendon grafts off muscle
- ▶ General orthopedics - periosteum elevator and spike

Available with or without side slot.

#### PRODUCT NO'S:


1114 [No Slot]  
Overall Length: 9.75" (24.8 cm)  
Handle Length: 4.25" (10.8 cm)  
Suture Hole: 2.5 mm x 13 mm

1114-01 [With Slot]  
Overall Length: 9.75" (24.8 cm)  
Handle Length: 4.25" (10.8 cm)  
Suture Hole: 2.5 mm x 13 mm



Advantage of the open slot -

- ▶ Convenient feeding and removal of sutures from slot
- ▶ Feeding of multiple thick sutures & sutures with knots
- ▶ Engaging and shuttling grafts with short suture loop ends



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362