



Kopplin Osteophyte Rongeur for Direct Anterior THA

Designed by Matthew Kopplin, MD

Designed to help remove osteophytes around the acetabulum in anterior THA

Helps to allow clean sharp cuts with better control. Thin flat tip helps to pass along the bone easier.



PRODUCT NO:

1771-04

Overall Length: 13" (35 cm)
Jaw Width: .375" (0.95 cm)
Jaw Length Bottom: .75" (1.9 cm)
Jaw Length Top: .625" (1.6 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



www.innomed.net

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.