

Lachiewicz Total Knee Revision Set

Designed by Paul F. Lachiewicz, MD

Used for total knee revision



PRODUCT NO.'S:

3700-00 [Complete Set In Case]

Individual Instruments:

Handle Length: 5.5" (14 cm)

3700-01 [10 mm Offset Edge Cutting Cement Chisel]

Chisel Offset: 10 mm

Chisel Length: 15 mm

Chisel Width: 8 mm

Length from Handle: 2.625" (6,7 cm)

3700-02 [15 mm Offset Edge Cutting Cement Chisel]

Chisel Offset: 15 mm

Chisel Length: 15 mm

Chisel Width: 8 mm

Length from Handle: 2.625" (6,7 cm)

3700-03 [Offset Femoral Component Dis-impactor]

Dis-impactor Width: 14 mm

Length from Handle: 3.375" (8,6 cm)

3700-04 [8 mm Cement Prosthesis Osteotome]

Osteotome Width: 8 mm

Blade Length: 2.25" (5,7 cm)

3700-05 [10 mm Cement Prosthesis Osteotome]

Osteotome Width: 10 mm

Blade Length: 2.25" (5,7 cm)

3700-06 [13 mm Cement Prosthesis Osteotome]

Osteotome Width: 13 mm

Blade Length: 2.25" (5,7 cm)

3700-07 [20 mm Cement Prosthesis Osteotome]

Osteotome Width: 20 mm

Blade Length: 2.25" (5,7 cm)

3700-08 [V-shaped Cement Splitter]

Splitter Width: .25" (.635 cm)

Blade Length: 2" (5,1 cm)

3700-09 [One-sided Cement Splitter]

Splitter End Overall Width: 9,3 mm

Splitter End Sharp Width: 6 mm

Splitter End Blunt Width: 3,3 mm

Blade Length: 3" (7,6 cm)

3700-10 [8 mm Cement Hook]

Hook Blade Width: 8 mm

Hook Shaft Diameter: 4 mm

Shaft Length: 5.375" (13,7 cm)

3700-11 [Cement Punch]

Punch Length: 3" (7,6 cm)

Punch Shaft Diameter: .25" (.63 cm)

3700-12 [Removal Cross Bar]

Overall Length: 4.375" (11,1 cm)

Shaft Diameter: 8 mm

3700-CASE [Case for Set]

Case Dimensions: 21" x 10" x 2.5" (53,3 x 25,4 x 6,35 cm)



USA MADE



10 mm Offset Edge Cutting Cement Chisel
#3700-01

15 mm Offset Edge Cutting Cement Chisel
#3700-02

Offset Femoral Component Disimpactor
#3700-03

8 mm Cement Osteotome
#3700-04

10 mm Cement Osteotome
#3700-05

13 mm Cement Osteotome
#3700-06

20 mm Cement Osteotome
#3700-07

V-shaped Cement Splitter
#3700-08

One-sided Cement Splitter
#3700-09

8 mm Cement Hook
#3700-10

Cement Punch
#3700-11

Removal Cross Bar
#3700-12

www.innomed.net



INNOMED

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2025 Innomed, Inc.