


Lambotte Osteotomes with Handle

Designed by John Cherf, MD

Designed with a handle for better control, which helps reduce rotation of the osteotome during use

The handle also provides a larger striking area for use with a mallet. The osteotome shafts are manufactured with stainless steel and are available both straight and curved in 1/4" blade width.



Straight #5260-01

Curved #5250-01

PRODUCT NO'S:

5260-01 [Straight]
Overall Length: 14" (35,6 cm)
Handle Length: 4.5" (11,4 cm)
Blade Width: .25" (6,3 mm)

5250-01 [Curved]
Overall Length: 14" (35,6 cm)
Handle Length: 4.5" (11,4 cm)
Blade Width: .25" (6,3 mm)

