


# Large Pin Distractor and Compressor

Larger 1/8" (3,2 mm) pin hole size for extra sturdy distraction or compression



LARGE PIN  
COMPRESSOR

LARGE PIN  
DISTRACTOR

**PRODUCT NO'S:**

4233 [Large Pin Distractor]  
Hole Diameters: For .125" (3,2 mm) K-wire Pins  
Overall Length: 8" (20,3 cm)

4234 [Large Pin Compressor]  
Hole Diameters: For .125" (3,2 mm) K-wire Pins  
Overall Length: 8" (20,3 cm)



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362