


# Lawrence Revision Knee Gap Balancing Tensioner Set

Designed by Jeffrey M. Lawrence, MD

*Designed to help tense the medial and lateral ligaments for gap balancing during revision surgery so that the AP cutting block does not impinge on the spreader during balancing*



<b>PRODUCT NO'S:</b>
1896-01 [Set - Left & Right]
<b>Also available individually:</b>
1896-01L [Left]
Overall Length: 9.25" (23,5 cm)
Pad Diameter: 1" (2,5 cm)
1896-01R [Right]
Overall Length: 9.25" (23,5 cm)
Pad Diameter: 1" (2,5 cm)



[www.innomed.net](http://www.innomed.net)



FREE TRIAL ON MOST INSTRUMENTS

**INNOMED**

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.