

Lee Buttress Pin Guides

Designed by Wonyong Lee, MD

Designed with pin guide exits acutely angled at 30° to help prevent the guide from slipping off the bone when placing the guide pin



PRODUCT NO'S:

3024-00 [Set of Three]

Set Includes/ Available Individually:

3024-01 [1.2/1.6 mm]

Overall Length: 6.125" (15.6 cm)

Handle Length: 4.875" (12.4 cm)

Guide Tube Length: 2" (5.2 cm)

Guide Angle from Handle: 30°

3024-02 [1.4/2.0 mm]

Overall Length: 6.125" (15.6 cm)

Handle Length: 4.875" (12.4 cm)

Guide Tube Length: 2" (5.2 cm)

Guide Angle from Handle: 30°





3024-03 [2.8/3.2 mm]

Overall Length: 6.125" (15.6 cm)

Handle Length: 4.875" (12.4 cm)

Guide Tube Length: 2" (5.2 cm)

Guide Angle from Handle: 30°



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362