


Levy Anterior Glenoid Retractor

Designed by Jonathan Levy, MD

Designed to help alleviate tension on anterior glenoid structures, while the handle is designed to optionally be clamped to the drape



PRODUCT NO:	
4536 [Small]	Overall Length: 10.5" (26,7 cm) Depth from Bend: 5.875" (14,9 cm) Blade Width: .75" (1,9 cm) Tooth Gap: .325" (8,2 mm)
4536-01 [Medium]	Overall Length: 10.5" (26,7 cm) Depth from Bend: 5.875" (14,9 cm) Blade Width: 1" (2,54 cm) Tooth Gap: .42" (10,6 mm)
4536-02 [Large]	Overall Length: 10.5" (26,7 cm) Depth from Bend: 5.875" (14,9 cm) Blade Width: 1.25" (3,2 cm) Tooth Gap: .55" (14 mm)



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.