



Levy Humeral Stem Extraction Punch

Designed by Jonathan Levy, MD

Ultra hard cobalt chrome shaft and impactor tip designed to help remove a humeral stem during revision total shoulder arthroplasty

Can be used to open up distal cement mantle or pedestal during revisions.



PRODUCT NO:

8627

Overall Length: 12" (30,5 cm)

Handle/Platform Length: 4.75" (12,1 cm)

Punch Rod Length: 7.25" (18,4 cm)

Platform: 3" x .75" (7,6 cm x 1,9 cm)

Shaft Diameter: 8 mm, tapers to 4 mm at tip



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362