


Light Wands – Short and Small Diameter

Short wand useful for proximal illumination, while thin diameter wands help illuminate deep cavities such as the femoral shaft

Can be attached to a fiber optic light cable with ACMI (female) connector.



Short Wand

Short Thin Wand

Long Thin Wand



Light Source Cable Adapters

PRODUCT NO'S:

8009-S [ACMI to Storz Adapter]
8009-W [ACMI to Wolf Adapter]



PRODUCT NO:

8011-01 [Short Light Wand]
Wand Length: 3" (7,6 cm)
Wand Shaft Diameter: 4.6 mm

8011-02 [Short Thin Light Wand]
Wand Length: 5.5" (14 cm)
and Shaft Diameter: 3.2 mm

8011-03 [Long Thin Light Wand]
Wand Length: 11" (28 cm)
Wand Shaft Diameter: 3.2 mm

