

Lighted Yankaur Suction Device

Designed by Adolph V. Lombardi Jr., MD

Designed to help provide effective suction with the addition of a light source for enhanced visualization

- ▶ Can be attached to a fiber optic light cable with ACMI (female) connector
- ▶ The handle is made of Delrin
- ▶ Entire device is steam sterilizable

PRODUCT NO:


8016-L

Overall Length: 11.75" (29,8 cm)

Handle Length: 3.93" (10 cm)

Handle Width: .86" (2,2 cm)

Suction Tube Diameter: .25" (6,35 mm)



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.