



## Standard Handle

# Lombardi Femoral Tibial Spreader

Designed by Adolph V. Lombardi Jr., MD

*Thin pads help to separate the femur and tibia during total knee procedures*



*The calibrated handle of the spreader helps to accurately gauge the gap, and makes it possible for two spreaders to be used to assist in balancing ligaments.*

### PRODUCT NO'S:

1875 [Large – Grooved Pads]

Overall Length: 9.25" (23,5 cm)

Pads: 22 mm x 13 mm

Opens to: 50 mm

1875-D [Large – Diamond Cut Pads]

Overall Length: 9.25" (23,5 cm)

Pads: 22 mm x 13 mm

Opens to: 50 mm

1876 [Small – Grooved Pads]

Overall Length: 7" (17,8 cm)

Pads: 22 mm x 13 mm


Opens to: 35 mm

1876-D [Small – Diamond Cut Pads]

Overall Length: 7" (17,8 cm)

Pads: 22 mm x 13 mm

Opens to: 35 mm



Grooved Pads Diamond Cut Pads



Large



Small



Large with Locking Mechanism

*Locking ratchet mechanism helps prevent accidental release, and provides for controlled adjustment and easy release.*



Small with Locking Mechanism

## Locking Mechanism

### PRODUCT NO'S:

1875-LR [Large – Grooved Pads with Locking Mechanism]

Overall Length: 9.625" (24,4 cm)

Pads: 22 mm x 13 mm


Opens to: 45 mm

1876-LR [Small – Grooved Pads with Locking Mechanism]

Overall Length: 7" (17,8 cm)

Pads: 22 mm x 13 mm

Opens to: 35 mm



## Small Grip Handle

### PRODUCT NO:

1876-SG [Small with Small Grip and Grooved Pads]

Overall Length: 7" (17,8 cm)

Pads: 22 mm x 13 mm

Opens to: 35 mm



*Designed with the grip closer together for easier gripping and to help reduce hand fatigue*



Small with Small Grip

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362

© 2022 Innomed, Inc.