





Lombardi Taper Cleaner

Designed by Adolph V. Lombardi Jr., MD

Designed to help clean a hip stem taper of corrosive byproducts prior to placement of the new femoral head



PRODUCT NO'S:	
Overall Length: 2.125" (5.4 cm)	
Outside Diameter: 1" (2.54 cm)	
8034	Small Short Taper 11.3/12.2 mm
8034-01	Long Taper 11.4/13.4 mm
8035-01	11/13 mm
8035-02	12/14 mm
8035-03	14/16 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362