



Lonner Swan Lateral Knee Retractor

Designed by Jess Lonner, MD and Martin Hyneman


Ergonomically designed for more effective retraction when using a robotic arm, allowing for clearer views of the surgical site

The retractor can effectively protect the lateral soft tissues and the patella when resecting the tibial plateau and lateral femoral condyle.

PRODUCT NO:

6651

Overall Length: 12.85" (32.6 cm)
 Blade Width at Body: 1" (2.54 cm)
 Blade Width Narrows To: .6" (1.5 cm)
 Blade Width at Tip: .13" (3.3 mm)



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
 info@innomed.net

912.236.0000 Phone
 912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
 info@innomed-europe.com

1.800.548.2362

© 2024 Innomed, Inc.