

Lotke Double Action Cartilage Graspers

Designed by Paul Lotke, MD

Double action strength helps to securely hold soft tissues


Angled to simulate the pinch forceps position. Ferris-Smith tips effectively hold soft tissues or needles. Powergrip avoids fatigue or excessive forces on the surgeons thumbs.



PRODUCT NO'S:

1710 [Standard]
Overall Length: 7.5" (19,1 cm)
1715 [Ratcheted]
Overall Length: 7.5" (19,1 cm)

MADE EXCLUSIVELY
FOR INNOVATED IN
GERMANY



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362