



Minimally Invasive Total Hip Surgery Retractors

Designed to help reduce soft tissue trauma during single incision minimally invasive total hip replacement surgery

Designed by
S. David Stulberg, MD




Narrow Cobra Retractors

Product No's:
6120 [Standard]
Overall Length: 11.75"
Blade Width: 16mm
6120-04 [Extra Long]
Overall Length: 15.5"
Blade Width: 16mm



Narrow Proximal Femoral Elevator

Product No's:
3420-08 [Narrow Prongs]
Overall Length: 11.75"
Blade Width at Widest: 45mm
3420-05 [Standard Prongs]
Overall Length: 11.5"
Blade Width at Widest: 45mm



Narrow Double Prong Acetabular Retractor

Product No: 6175
Overall Length: 12.5"
Blade Width: 32mm



MIS Hip Retractors

Product No's:
6260 [Extra Narrow]
Overall Length: 12"
Blade Width: 20mm
Handle Length: 9.5"
Depth from Bend: 6"
6265 [Extra Long Extra Narrow]
Overall Length: 15"
Blade Width: 20mm
Handle Length: 11"
Depth from Bend: 8"



Narrow Assistant Free™ Blades

Product No's:
7480-02 2" Blade Depth
7480-03 3" Blade Depth



Single Prong Double Bent Hohmann Acetabular Retractor – Extra Long

Product No's:
6214 [Longer Blade]
Overall Length: 18"
Blade and Tip Length: 5"
Blade Width: 15mm
6210-04 [Standard Blade]
Overall Length: 16.25"
Blade and Tip Length: 3"
Blade Width: 15mm






Single & Double Prong Double Bent Hohmann Acetabular Retractor – Long

Product No's:
6213 [Single Prong – Longer Blade]
Overall Length: 15"
Blade and Tip Length: 5"
Blade Width: 15mm
6210-02 [Single Prong – Standard Blade]
Overall Length: 12.5"
Blade and Tip Length: 3"
Blade Width: 15mm
6220 [Double Prong – Standard Blade]
Overall Length: 12.5"
Blade and Tip Length: 3"
Blade Width: 15mm



Single Prong Double Bent Hohmann Acetabular Retractor

Product No's:
6212 [Longer Blade]
Overall Length: 11.25"
Blade and Tip Length: 3.5"
Blade Width: 15mm
6210 [Standard Blade]
Overall Length: 10.5"
Blade and Tip Length: 2.5"
Blade Width: 15mm



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone **800 548 2362**
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

