



Helps provide excellent visibility and ligament protection during Total and Unicompartmental Knee Replacement Surgery

Minimally Invasive Knee Retractors


PRODUCT NO'S:	
1	S3035 [Small Hohmann Retractor] Overall Length: 7.5" (19,1cm) Blade Width: 25mm
2	S3036 [Large Hohmann Retractor] Overall Length: 8" (20,3cm) Blade Width: 36mm
3	S3037 [Condylar Retractor] Overall Length: 7.5" (19,1cm) Blade Width: 12mm
4	S3038 [Superior Retractor] Overall Length: 8.75" (22,2cm) Blade Width: 31mm
5	S3039 [Patellar Retractor] Overall Length: 10.25" (26cm) Blade Width: 45mm
6	S3042 [Soft Tissue Retractor] Overall Length: 8.75" (22,2cm) Blade Width: 36mm



1 MI Small Hohmann Retractor



4 MI Superior Retractor




2 MI Large Hohmann Retractor



5 MI Patellar Retractor



3 MI Condylar Retractor



6 MI Soft Tissue Retractor

Retractor handles are designed for optional use with modular weights to help hold the retractors in place.



Modular Weights

PRODUCT NO'S:	
3430-01	1.5 lbs. (.68 kg)
3430-02	2.0 lbs. (.91 kg)
3430-03	2.5 lbs. (1.13 kg) with attaching hook

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to Launch Our Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362