



Minimally Invasive Knee Retractors

Helps provide excellent visibility and ligament protection during Total and Unicondylar Knee Replacement Surgery

PRODUCT NO'S:

S3035 [Small Hohmann Retractor]

Overall Length: 7.5" (19,1cm)

Blade Width: 25mm

S3036 [Large Hohmann Retractor]

Overall Length: 8" (20,3cm)

Blade Width: 36mm

S3037 [Condylar Retractor]

Overall Length: 7.5" (19,1cm)

Blade Width: 12mm

S3038 [Superior Retractor]

Overall Length: 8.75" (22,2cm)

Blade Width: 31mm

S3039 [Patellar Retractor]

Overall Length: 10.25" (26cm)


Blade Width: 45mm



1 MI Small Hohmann Retractor



4 MI Superior Retractor



2 MI Large Hohmann Retractor



5 MI Patellar Retractor



3 MI Condylar Retractor

Retractor handles are designed for optional use with modular weights to help hold the retractors in place.



Modular Weights

PRODUCT NO'S:

3430-01 1.5 lbs. (.68 kg)

3430-02 2.0 lbs. (.91 kg)

3430-03 2.5 lbs. (1.13 kg) with attaching hook

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362