




# Macko Square Tipped Rongeur

Designed by Victor W. Macko, MD

*Square tipped rongeur with an ergonomic grip, double action mechanism, long reach, and low profile – for use in total ankle, hip, knee, and spine surgery*

In total ankle surgery, the aggressive, low profile jaws aid in the removal of tibial bone in spite of limited space. The square ended tip helps produce a flat, finished surface following anterior talar facet reaming.

When used for morcelizing bone graft, the shallow, wide jaw helps avoid impaction.



Three Jaw Sizes Available

PRODUCT NO'S:	
1778-01	Jaw Bite: 5 x 18mm Overall Length: 10"
1778-02	Jaw Bite: 7 x 18mm Overall Length: 10"
1778-03	Jaw Bite: 10 x 18mm Overall Length: 10"

**PROUDLY  
MADE IN THE USA**



ISO 9001:2008 • ISO 13485:2003

**FREE TRIAL ON MOST INSTRUMENTS**

Scan to  
Launch Our  
Website



**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

**1.800.548.2362**

© 2013 Innomed, Inc.