



Mazzara Rongeur for Small Bones

Designed by James T. Mazzara, MD

Designed for bone and soft tissue removal in small joint surgery, the pistol grip handle lessens hand fatigue and slippage, and allows for better visualization
Available in two jaw sizes.



PRODUCT NO'S:	
1765-04	[2 x 10 mm Jaw Bite] Overall Length: 7.25" (18,4 cm)
1765-05	[4 x 10 mm Jaw Bite] Overall Length: 7.25" (18,4 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362