

McFarland Bent Cobb Elevator

Designed by Edward McFarland, MD

Designed for retraction while helping to protect the axillary nerve in shoulder surgery

Ultra hard titanium nitride coating helps to help to prolong sharpness.



PRODUCT NO:

3431

Overall Length: 9.5" (24.1 cm)

Length from Bend: 3.5" (8.9 cm)

Cobb End Width: .8" (20 mm)

Angle of Bend: 30°



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362