


McMaster Abductor Retractor

Designed by William D. McMaster, MD

Designed to help with proximal femur exposure helping to protect the abductors – gluteus medius and minimus – during posterior approach THA

The ergonomic design allows application where soft tissue retraction is needed.



PRODUCT NO:
6385
Overall Length: 8.75" (22,2 cm)
Handle Length: 4.625" (11,7 cm)
Blade Depth: 4.125" (10,8 cm)
Blade Width: 1.5" (3,8 cm)
Blade Curve Diameter: 1.75" (4,4 cm)



www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362