

Medial Acetabular Retractors with Large Handle

Designed for acetabular exposure during total hip surgery



PRODUCT NO'S:

C1001 [Left]

Overall Length: 15" (38,1 cm)

Handle Length: 6" (15,2 cm)

Blade Width: 32 mm

C1002 [Right]

Overall Length: 15" (38,1 cm)

Handle Length: 6" (15,2 cm)

Blade Width: 32 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED