



# Medial Malleolus Fracture Reduction Aid

Designed by Christopher Blair, DO

*Designed to hook under the medial malleolus to help reduce the medial malleolus fragment while two K-wire guides supply trajectory for wires*  
For K-wires up to 1,6 mm (.062")



**PRODUCT NO:**  
3664  
Overall Length: 7" (17,8 cm)  
Handle Length: 4" (10,2 cm)  
Plate Width: .8" (2 cm)  
Plate Length: 3" (7,6 cm)  
Guide Tube Length: 6 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe

Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362

# INNOMED