

Mehalik Posterior Glenoid Retractor with Long Handle

Designed in collaboration with Mayo Clinic, modified by John Mehalik, MD.

Designed to help expose the posterior aspect of the glenoid



PRODUCT NO:

1909
Overall Length: 13.2" (33,5 cm)
Access Hole Internal Diameter: 36 X 30 mm



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362