

Modified Mini Hohmann Retractors

Designed by Jeffrey Lawton, MD

Used for small bone surgery



PRODUCT NO'S:	
1665	[Blade: 6 mm Wide / 35 mm Drop] Overall Length: 5.875" (14,9 cm) Blade Width: 6 mm Blade Drop: 35 mm
1665-01	[Blade: 6 mm Wide / 17 mm Drop] Overall Length: 5.5" (14 cm) Blade Width: 6 mm Blade Drop: 17 mm
1666	[Blade: 8 mm Wide / 35 mm Drop] Overall Length: 5.875" (14,9 cm) Blade Width: 8 mm Blade Drop: 35 mm
1666-01	[Blade: 8 mm Wide / 17 mm Drop] Overall Length: 5.5" (14 cm) Blade Width: 8 mm Blade Drop: 17 mm
1666-02	[Blade: 8 mm Wide / 17 mm Drop] Overall Length: 6.25" (15,9 cm) Blade Width: 8 mm Blade Drop: 17 mm
1666-LG	[Blade: 8 mm Wide / 72 mm Drop] Overall Length: 7.125" (18,1 cm) Blade Width: 8 mm Blade Drop: 72 mm

New!
New!

New!

New!

www.innomed.net

FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.