


Mobasser Bayonet Forceps

Designed by Jean Pierre Mobasser, MD

Tip design can hold a larger amount of bone graft, and then be used as a tamp to push it down into the disc space

Bayonetted forceps has a long reach, which is particularly helpful in less invasive spine surgery.



PRODUCT NO:

5082

Overall Length: 10"

Forceps Ends: 5x10mm

**PROUDLY
MADE IN
THE USA**

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED 

1.800.548.2362