



Modified Angled Hohmann Retractor with Long Handle and Short Tip

Designed by R. Michael Meneghini, MD

Longer handle to help provide safe patella retraction with excellent ergonomics, and useful in other orthopedic procedures

Excellent for gently retracting the patella and extensor mechanism on heavy obese patients and muscular male patients.



PRODUCT NO:
7119
Overall Length: 11.5" (29.2 cm)
Handle Length: 9.25" (23.5 cm)
Blade Depth: 5" (12.7 cm)
Blade Width: .75" (19 mm)
Prong Depth: 1" (2.54 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362