

Modified Darrach-type Bent Elevator

Design modification by R.L. Stowell, MD of original design by Evan Flatow, MD

Designed for difficult glenoid exposure, the elevator is placed around the posterior glenoid rim, retracting the cut humeral surface

The retractor is serrated on both sides to help prevent slippage.

PRODUCT NO:

1966

Overall Length: 10" (25,4 cm)

Blade Depth: 5" (12,7 cm)

Blade Width: 1" (2,54 cm)



◀ Serrated on both sides

www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362