


# Modified Wide Hohmann Retractor with Taylor Tip


Designed by Jeffrey P. Beckenbaugh, DO

*Anterior and posterior acetabular retractors for all approaches, including the direct anterior approach, featuring a hammer platform for insertion with a mallet*

Used as a calcar and posterior femoral retractor for the posterior approach, and an anterior femoral elevator for the direct anterior approach.



|                                     |
|-------------------------------------|
| <b>PRODUCT NO'S:</b>                |
| 3012-00 [Set - Left & Right]        |
| <b>Also available individually:</b> |
| 3012-L [Left]                       |
| Overall Length: 14.75" (37,5 cm)    |
| Blade Depth: 5" (12,7 cm)           |
| Blade Width: 45 mm                  |
| 3012-R [Right]                      |
| Overall Length: 14.75" (37,5 cm)    |
| Blade Depth: 5" (12,7 cm)           |
| Blade Width: 45 mm                  |



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362

© 2022 Innomed, Inc.