




Modified Humeral Head Retractors

Designed to help lever and displace the proximal humerus posteriorly

PRODUCT NO:

2436 [Deep]
Overall Length: 10" (25,4 cm)
Blade Width: 26 mm

2437 [Shallow]
Overall Length: 11" (30 cm)
Blade Width: 26 mm



www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362