



Modified Lambotte Cup Removal Osteotomes

Designed with different hemisphere of curves to match cups of different sizes

Four osteotomes with different hemispherical radii allow the osteotomes to fit next to the outer surface of different size acetabular hip cups. The osteotomes have a handle for better control and a hammering platform.



PRODUCT NO'S:	
5240-44	Blade Radius: 44 mm Overall Length: 12.75" (32,4 cm) Handle Length: 4.75" (12 cm)
5240-48	Blade Radius: 48 mm Overall Length: 12.75" (32,4 cm) Handle Length: 4.75" (12 cm)
5240-52	Blade Radius: 52 mm Overall Length: 12.75" (32,4 cm) Handle Length: 4.75" (12 cm)
5240-56	Blade Radius: 56 mm Overall Length: 12.75" (32,4 cm) Handle Length: 4.75" (12 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.