




# Modified Smith-Peterson Style Osteotomes for Acetabular Cup Removal

Designed by Merrill Ritter, MD

*Multi-arch osteotomes help in removal of total hip cups*

For removal of total hip cups, the different curvatures help to fit next to a cups outer surface. The osteotomes have a handle for better control, plus a hammering platform end.



#### PRODUCT NO'S:

5280-01 [Short]

Blade Dimensions: 20 mm x 20 mm

Overall Length: 10.875" (27,6 cm)

Handle Length: 5" (12,7 cm)

5280-02 [Medium]

Blade Dimensions: 20 mm x 35 mm

Overall Length: 11.675" (29,6 cm)

Handle Length: 5" (12,7 cm)

5280-03 [Long]

Blade Dimensions: 20 mm x 50 mm

Overall Length: 12.25" (31,1 cm)

Handle Length: 5" (12,7 cm)



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.